

STROBOSCOPE

STROBOLIGHT

EXTRA BRIGHT

350W XENON LAMP

LONG-LIFE LAMP MAKES THE STROBOLIGHT VIRTUALLY MAINTENANCE FREE



ECLERIS  TECHNOLOGY

- AUTOMATIC VOICE SYNCHRONIZATION
- TOTAL FOOTSWITCH CONTROL



ECLERIS

SUITABLE FOR BOTH FLEXIBLE
AND RIGID ENDOSCOPES



TECHNICAL SPECIFICATIONS

MODEL	STROBOLIGHT
WORKING RANGE	78 - 1100 HZ
DISPLAY RANGE	120 - 999 HZ
AUTO MODE	Large / Fine Tune. Full automatic fundamental frequency detection controlled by double footswitch pedal for speed control and phase position.
MANUAL MODE	Manual fundamental frequency detection controlled on frontal panel display.
OPERATING MODES	Slow Motion Mode Freeze Frame Mode
MICROPHONE	Direct Contact
WATTAGE	350 Watts Xenon
LAMP LIFE	> 2000 Hours
STROBOSCOPIC LIGHT	Xenon
PILOT LIGHT	Xenon for no color change
VOLTAGE RANGES (SELECTABLE)	230 V ± 10 %, 50 / 60 Hz 115 V ± 10 %, 50 / 60 Hz
MAIN UNIT	15,35" x 11,81" x 5,7" 39 x 30 x 14.5 cm
PEDAL	11,02 x 10,82" 28 x 27.5 cm
MAIN UNIT	17,63 Pounds 8 Kgs
PEDAL	8,81 Pounds 4 Kgs

4.05 / 03 / 2007

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE.



DISTRIBUTED BY:



ECLERIS
www.ecleris.com