

PICTURE / AUDIO / VIDEO



- PORTABLE
- USER FRIENDLY



- PORTABLE
- USER FRIENDLY
- EXTERNAL USB
- FOOTSWITCH
- INPUT / OUTPUT
- PATIENT DATA BASE
- PATIENT REPORTS WITH PICTURES
- AUDIO AND VIDEO FILES



MINIMUM REQUIREMENTS FOR RUNNING ENDODIGI SOFTWARE

OPERATING SYSTEM: - WINDOWS XP OR WINDOWS VISTA
 PROCESSOR: - PENTIUM IV 2.2 GHZ OR FASTER
 - ATHALON 64 2.0 GHZ OR FASTER
 HARD DRIVE: - 80Ghz (or LARGER)
 PORTS USB 2.0: - 2 AVAILABLE
 SCREEN RESOLUTION: - 1024 X 768
 CD-RW: - RECOMMENDED FOR BACKUP

TECHNICAL SPECIFICATIONS

| | |
|--------------|---|
| MODEL | ENDODIGI |
| CLASS | 1 |
| WEIGHT | 1.5 Kg 3 lb. |
| DIMENSIONS | 12 x 4 x 25 cm. 4.72 x 1.57 x 9.84" |
| NORMS | PAL NTSC |
| RESOLUTION | 24 BITS (RGB) |
| FILE FORMATS | JPEG (PICTURES) WMV9 MPEG4 ENCODER (VIDEO) |
| INTERFACE | USB 2.0 |

DISTRIBUTED BY:

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE. FL-20001 rev.10.06.08