

NASOPHARYNGOSCOPES

The new generation of nasopharyngoscopes convince by their intuitive operability and optimal ergonomics, designed for the tactile sensation and non-tiring handling as required in daily routine operation. High-End fibre-optic bundles and a high capacity system produce large and bright images at best resolution – necessary for save diagnosis and therapy.

Due to the consequent use and reliable processing of high-class materials, the nasopharyngoscopes are robust and have a long lifetime in spite of the increasing stress of modern cleaning and sterilizing facilities.

Artikel-Nr. >



Spezifikationen

- > functional design for easy and comfortable handling
- > Sensitive yet robust control
- > Optimal flexibility and torsional stiffness of the insertion channel
- > high-performance optical system for larger, brighter and sharper images
- > suitable reprocessing procedures:
 - fully immersible in cleaning and disinfectant solutions
 - STERIS®
 - ETO
 - machine reprocessing (RDG-E)

ACUTRONIC Medical Systems AG

Fabrik im Schiffli | 8816 Hirzel | Switzerland | Phone ++41 44 729 70 80 | Fax ++41 44 729 70 81
info@acutronic-medical.ch | www.acutronic-medical.ch



	EF-N 130 400 034	EF-NS 130 400 134	EF-N slim 130 400 028	EF-N14 130 414 040	EF-N14 S 130 414 140
diameter of working channel	--	--	--	1.4 mm	1.4 mm
incorporated light-guide cable	--	x	--	--	x
sheath diameter	3.4 mm	3.4 mm	2.8 mm	4.0 mm	4.0 mm
working length	320 mm	320 mm	320 mm	320 mm	320 mm
bending up/down	130°	130°	130°	130°	130°
bending radius	8 mm	8 mm	8 mm	8 mm	8 mm
optical system					
field of view	80°	80°	80°	80°	80°
depth of field	1-50 mm	1-50 mm	1-50 mm	1-50 mm	1-50 mm