

PROCAM SD

VIDEO CAMERA

- ENDOSCOPY / COLPOSCOPY / MICROSCOPY



- High resolution 500 lines
- All controls on handpiece
- White balance: Anti-moiré / Windows / Photo
- Universal CCD 1/3"

PROCAM SD

VIDEO CAMERA



POWER SUPPLY




NTSC 525 LINES

OPTIONAL LIGHT SOURCES

EVERLUX PRO



XENOLUX180



TECHNICAL SPECIFICATIONS

Resolution	Pal 625 lines / ntsc 525 lines
Ccd	1/3"
Shutter speed	Up to 1/100.000 sec.
Min. Illumination	2 lux f1.4
Video output	Composite / S-VIDEO
Power supply	100v / 240v
Frequency	50 - 60hz

Lamp	LED
Light temp. Color	6000 K
Fiber optic connection	Touret compatible with main brands
Power required	110 w
Power supply	100v / 240vca
Frequency	50 - 60hz

Lamp	Xenon 180 w
Light temp. Color	6000° kelvin
Fiber optic connection	Touret compatible with main brands
Power required	250 w
Power supply	100v / 240v
Frequency	50 - 60hz

ECLERIS RESERVES THE RIGHT TO MODIFY THIS PRODUCT WITHOUT PRIOR NOTICE

Rev. 22.12.14



www.ecleris.com

ECLERIS
INNOVATION VALUE