

RSX **USB**

Video Rhino Laryngoscope

Ultra-compact with integrated LED and USB connection



First ultra-compact scope in the world | Integrated high power LED
Advanced Image Technology (AIT) | Best price performance ratio
Reduced weight and improved mobility | USB - one cable to use all functions



Chip on the tip
with precise optics

integrated
high power LED

interface

power on/off



angulation

capture button

power LED



leakage test



Technical Data

RSX USB Video Rhino Laryngoscope

Diameter	3.9 mm
Working length	310 mm
Tip Articulation	130°/130°
FOV	90°
Resolution	400 x 400 pixel / 160.000 pixel (PAL)
Video interface	USB
Illumination	integrated High Power LED
Protection class (environmental influences)	IP68
Protection class (electric shock)	class II
Risk class (German Medical Devices Act)	1