

Naso-Pharyngoscope



Flexible precision-endoscope with high performance optic.

Engineered for highest standards and daily requirements. A long life period is guaranteed by the high resolved bundle-light, combined with mechanical components of high precision and a sturdy construction. The highest quality requirements are complied by the bundle-light with 10.000 pixels and the multi-lenses objective.



Technical data

Image bundle 10.000 pixel Min.

Working length 300 mm

Working Diameter 3,8 mm

Flexibility of distal end 2 x 125°

Length of bendable distal end 25 mm

Min. bending-radius 30 mm

Viewing direction 0°

Viewing angle 90°

Depth of focus 5 mm until infinity

Weight 255 g

Adapter types (included): ACMI, Wolf and Storz

Disinfection/Sterilization

In liquid disinfectant
Gas-sterilizable

Order numbers

Naso-Pharyngoscope 3,8 mm incl. case
Test device to verify sealing of flexible tube

no. 3900
no. 1360.007

Rudolf Riester GmbH & Co. KG
P.O. Box 35 • DE-72417 Jungingen
Germany
Tel.: +49 (0)74 77/92 70-0
Fax: +49 (0)74 77/92 70-70
info@riester.de
www.riester.de

Subject to alterations

ENT- Endoscopy