



LiteOptics Hysteroscopes

Specially designed for the ambulatory setting

Diagnostic Hysteroscope



Operating Hysteroscope



LiteOptics Ambulatory Hysteroscope

Please specify which items you require to make up your own personal sets.

1.9mm 0° SemiFlex HD Telescope 312mm

Code; URO.020.00.3120



1.9mm 25° SemiFlex HD Telescope 312mm

Code; URO.020.09.3120



2.9mm Single Flow Diagnostic Hysteroscope Sheath with Rotating Tap

Long working length 275mm

Code; 200-131



3.5mm x 4.2mm One Piece Continuous Flow Hysteroscope Sheath with 5 Fr Operating Channel

Long working length 235mm

Code; 200-373



3.5mm x 2.3m Light Guide Cable for High Intensity Light Sources

Code; 00-22330-40



Serrated Cup Grasping/Biopsy Forceps

Code; 56425S



Single Action Scissors


Code; 57508S



Sterilisation Tray with Lid

Suitable for one complete set comprising:
Telescope, diagnostic and operating sheath,
forceps, scissors and light cable

Code; CIL2701



VTW

Bipolar 5fr electrode

Equivalent to Twizzel electrode (Gynaecare ref: 00467)

Compatible with Versapoint 1 and 2



Versascope Adaptor

Converts LiteOptics 1.9mm telescopes to fit
Gynaecare Versascope sheath

Code; 200-378



Full repair service available

Certification

ISO 9001

EN ISO – 13485

CE 0483 EC 93/42/EEC

NHS Supply Chain approved



Tel; 01799542716

Email; enquiries@liteoptics.com